

PGxPredict Test

Sample Collection Instructions

Customer Service (Toll Free) 877.2.PGX.Health (877.274.9432)

For the healthcare provider:

The physician and patient should review, complete, and SIGN the PGxPredict Tests Requisition/ Informed Consent. Return the original copy to PGxHealth with the sample and retain a copy for the patient's medical record.



1
Write the patient's initials and date of birth on the sample collection cards. Be careful not to touch the sample collection area.

Complete the information requested on the outside of the white envelope, including the date of sample collection.



2
Clean a puncture site with the alcohol wipe. Have the sample collection cards ready.



3
Pull the yellow indicator tab off of the safety lancet. Position safety lancet firmly against finger. Hold lancet between fingers and place thumb on white button.

Press white button firmly into blue housing. Button will lock into housing and will not recoil. Do not pull lancet away from puncture site until after activation.



4
Gently massage from base of finger to near the puncture site, holding patient's hand below elbow level.



5
Please completely fill at least 2 sample collection spots on the card(s). Allow samples to dry for 20 minutes. Place in white envelope and seal.



6
Place the white envelope into the Infecon bag (verify that the desiccant is in the bag). Place the Test Requisition/Informed Consent in the outer sleeve of the Infecon bag.

Ensure all the information requested on the outside of the white envelope is complete.



7
Affix the preprinted label to the top of the shipping container and ship on the day of collection to:

PGxHealth
Attn: PGxPredict Tests
Five Science Park
New Haven, CT 06511